



**Lámparas / Lamps**

***L33-38 series***

***Disamed***

DISTRIBUIDORA ANDALUZA DE MEDICINA, S.L.  
Jaime Balmes, 6- 41007 Sevilla  
Tfnos.: 954574032 (5 líneas) Fax.: 954584979  
Web: [www.disamed.es](http://www.disamed.es) - correo: [info@disamed.es](mailto:info@disamed.es)



# Lámpara de examen modelo L33-38

## Examination lamp model L33-38

Lámpara indicada para cirugía menor y exploración avanzada. Sus óptimas dimensiones que permiten un manejo fácil y fiable, con una iluminación uniforme en área de trabajo, así como un alto índice de seguridad tanto para el paciente como para el usuario.

Alta rentabilidad por su simple mantenimiento, mínimo consumo y larga vida útil de sus componentes.

*Lamp suitable for minor surgery and advanced exploration. The optimal dimensions that allow easy and reliable to use, with a uniform illumination of the work field, and a high level of safety for the patient and the user.*

*High profitability for the simple maintenance, low consumption and long life of the components.*

DATOS TÉCNICOS		TECHNICAL DATA
Intensidad lumínica	100 cm : 70.000 lux	Illumination intensity
Temperatura de color	4.750 °K	Colour temperature
Ø campo de trabajo d10	35cm	Ø Light field d10
Ø campo de trabajo d50	18cm	Ø Light field d50
Profundidad campo (L1+L2)	>150cm	Depth of illumination
Módulos LED	8	Modul of LEDs
Vida media LED	40.000 h	Lifetime LED
Potencia	60 W	Power
Voltaje	100 – 240 V / 47-63 Hz	Voltage
Mango esterilizable	Si / Yes	Sterilizable handle
Tamaño de la cúpula	Ø 330 mm	Size of light head
Color	RAL 9002	Colour
Peso bruto / neto	L21-25R 19,5 Kg / 15 Kg	Gross Weight / Net
Dimensión embalaje	L21-25R 108 x 36 x 44 cm /0,171 m <sup>3</sup>	Packing size

MODELOS / MODELS  
L33-38 series

L33-38R	rodable <i>moveable lamp</i>
L33-38RE	rodable con equipo emergencia <i>moveable with emergency equipment</i>
L33-38P	con soporte a pared <i>wall mounted</i>
L33-38T	de techo <i>ceiling</i>



Tecnología, 2 – nave 7  
08780 – PALLEJA ( Barcelona ) – SPAIN  
Teléfono: (+34) 93 334 01 12  
Fax: (+34) 93 440 25 64  
<http://www.ordisi.com>  
[comercial@ordisi.com](mailto:comercial@ordisi.com)  
[export@ordisi.com](mailto:export@ordisi.com)

